


# Pt. III NMR SERVICE SAMPLE SUBMISSION FORM

Name	Supervisor & Group																																
Personal Code	Date																																
Name of compound	Lab																																
	Phone																																
	Email																																
<b>Sample Information</b>																																	
Structure/Molecular Weight	<b>Sample will not be run if the information below is not provided, or is in an incorrect tube.</b>		TUBE LENGTH MINIMUM 175mm																														
	<div style="display: flex; justify-content: space-between;"> <div style="width: 40%;">Solvent</div> <div style="width: 60%;"></div> </div>																																
<div style="display: flex; justify-content: space-between;"> <div style="width: 40%;">Sample Weight</div> <div style="width: 60%;"></div> </div>																																	
<table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th style="width: 20%;">1D EXP.</th> <th style="width: 20%;"></th> <th style="width: 20%;">2D EXP.</th> <th style="width: 20%;"></th> <th style="width: 40%;">Supervisor's Signature</th> </tr> </thead> <tbody> <tr> <td><sup>1</sup>H</td> <td></td> <td>COSY</td> <td></td> <td></td> </tr> <tr> <td><sup>13</sup>C</td> <td></td> <td>HSQC</td> <td></td> <td></td> </tr> <tr> <td>APT/DEPT</td> <td></td> <td>HMBC</td> <td></td> <td></td> </tr> <tr> <td></td> <td></td> <td>NOESY</td> <td></td> <td></td> </tr> <tr> <td></td> <td></td> <td></td> <td></td> <td></td> </tr> </tbody> </table>				1D EXP.		2D EXP.		Supervisor's Signature	<sup>1</sup> H		COSY			<sup>13</sup> C		HSQC			APT/DEPT		HMBC					NOESY							
1D EXP.		2D EXP.		Supervisor's Signature																													
<sup>1</sup> H		COSY																															
<sup>13</sup> C		HSQC																															
APT/DEPT		HMBC																															
		NOESY																															

## Safety and Hazard Information

Please indicate where appropriate

Carcinogenic	
Toxic Vapours	
Lachrymatory	

Poisonous	
Corrosive	
Stench	

Other Hazard

## Sample Submission Guide

Ensure your NMR tube is a Wilmad 528 tube of the correct length; non-standard tubes will not be run. Ensure you have filled in the safety information. Fill in the molecular weight and weight of sample so we can choose which instrument is best for you. Good results start with careful preparation – use the solvent volume depth guide above and avoid floating impurities. If you have a question ask one of the NMR Service team.